

Human PD-L1 (B7-H1) protein (Recombinant) {Biotin} (STJP000218)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant human PD-L1 (B7-H1) is a protein consisting of 225 amino acid residue subunits, due to glycosylation migrates as an approximately 40kDa protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
UniProt ID	PD1L1_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--