

Mouse Interleukin-10 protein (Recombinant) (STJP000235)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO Cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.0.
Molecular Weight	Native Mouse Interleukin-10 is generated by the proteolytic removal of the signal peptide and propeptide, the molecule has a calculated molecular mass of approximately 19kDa. Recombinant Mouse IL-10 is a protein consisting of 161 amino acid residue su
Storage Instruction	Can be stored in working aliquots at 4°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	16153
Gene Symbol	IL10
UniProt ID	IL10_MOUSE
Immunogen Region	Mature chain
Immunogen Sequence	MPGSALLCCL LLLTGMRISR GQYSREDNNC THFPVGQSHM LLELRTAFSQ VKTFFQTKDQ LDNILLTDSL MQDFKGYLGC QALSEMIQFY LVEVMPQAEK HGPEIKEHLN SLGEKLTLR MRLRRCHRFL KCENKSKAVE QVKSDFNKLE DGGVYKAMNE FDIFINCIEA YMMIKMKS NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--