

MERS Spike protein S2 protein (Recombinant) (STJP000255)

GENERAL INFORMATION

Product Type Proteins

Host / Source CHO Cells

PRODUCT PROPERTIES

Formulation Lyophilised from 0.2 μ m filtered PBS solution pH7.4.

Molecular Weight Recombinant MERS-CoV Spike S2 is a protein consisting of amino acid residues [748-1296], due to glycosylation migrates as an approximately 67 kDa protein on SDS-PAGE. NA

Storage Instruction The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Immunogen Region ECD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.