

2019 nCoV/COVID19 Nucleocapsid protein (Recombinant) (His-Tag) (STJP000259)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO Cells

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS, pH 7.4.
Molecular Weight	Recombinant 2019-nCoV Nucleocapsid protein is a protein consisting of 431 amino acid residues, due to glycosylation migrates as an approximately 70 kDa protein on SDS-PAGE. NA
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Immunogen Region	1-419 aa
------------------	----------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



1 μ g loaded Tris-Bis PAGE under reducing condition. The purity is greater than 95%.