

2019 nCoV/COVID19 Spike protein S protein (Recombinant) (STJP000265)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered 20 mM Tris, 300 mM NaCl, 10 % glycerol, pH 8.0.
Molecular Weight	Recombinant 2019-nCoV Spike S protein is a protein consisting of the S1+S2 subunits contains 1210 amino acid residues, due to glycosylation migrates as an approximately 120 kDa protein on SDS-PAGE. NA
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Immunogen Region	ECD
------------------	-----

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



1 μ g loaded Tris-Bis PAGE under reducing condition. The purity is greater than 95%.