

2019 nCoV/COVID19 Spike protein S1 protein (Recombinant) (His-Tag) (STJP000269)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant NCOV S1 encodes the NCOV spike glycoprotein S1 sequence. This molecule has a calculated mass of approximately 76kDa. Recombinant S1 migrates due to glycosylation as an approximately 120 kDa protein under reducing conditions in SDS-PAGE. N
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Immunogen Region	ECD
------------------	-----

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



1 μ g loaded Tris-Bis PAGE under reducing condition. The purity is greater than 95%.