

Human Angiopoietin-2 protein (Recombinant) (STJP000283)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Host / Source | HEK293 |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose. |
| Molecular Weight | Recombinant Angpt2is a protein consisting of490 amino acid residue subunits, due to glycosylation migrates as an approximately 70kDa protein on SDS-PAGE. NA |
| Storage Instruction | Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 285 |
| Gene Symbol | ANGPT2 |
| UniProt ID | ANGP2_HUMAN |
| Immunogen Region | ECD |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|