

Human ADAM 12 protein (Recombinant) (STJP000288)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant human ADAM12 (meltrin alpha) is a protein consisting of 490 amino acid residue subunits, due to glycosylation migrates as an approximately as 65kDa (proprotein) , 23kDa (propeptide) and 46kDa (mature form) band protein on reduced SDS-PAGE
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	8038
Gene Symbol	ADAM12
UniProt ID	ADA12_HUMAN
Immunogen Region	ARG29-ASP513

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--