

Human KLF4 protein (Recombinant) (STJP000289)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution pH 7.4.
Molecular Weight	Recombinant Human KLF4 is a monomer consisting of 470 amino acid residue subunits and migrates as an approximately 51 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	9314
Gene Symbol	KLF4
UniProt ID	KLF4_HUMAN
Immunogen Region	Mature chain

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--