

Human Interleukin-18 Binding Protein (IL-18BP) protein (Recombinant) (STJP000301)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	0.2 Mu m filtered PBS solution.
Molecular Weight	Native human IL-18BP, generated by the proteolytic removal of the signal peptide and propeptide, the molecule has a calculated molecular mass of approximately 19kDa. Recombinant Interleukin-18 Binding Protein is a monomer protein consisting of 163 am
Storage Instruction	Can be stored in working aliquots at°C-°C C for one month, or at-20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	10068
Gene Symbol	IL18BP
UniProt ID	I18BP_HUMAN
Immunogen Region	Mature chain
Immunogen Sequence	MTMRHNWTPD LSPLWVLLC AHVVTLLVRA TPVSQTTTAA TASVRSTKDP CPSQPPVFPA AKQCPALEVT WPEVEVPLNG TLSLSCVACS RFPNFSILYW LGNGSFIEHL PGRLWEGSTS RERGSTGTQL CKALVLEQLT PALHSTNFSCLVLDPEQVVQ RHVVLAQLWA GLRATLPPTQ EALPSSHSSP QQQG NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--