

Human Renin protein (Recombinant) (STJP000303)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Molecular Weight	Recombinant human Renin is a protein consisting of 390 amino acid residue subunits, due to glycosylation migrates as an approximately 45kDa protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	22943
Gene Symbol	DKK1
UniProt ID	DKK1_HUMAN
Immunogen Region	LEU24-ARG406, ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--