

Human TGF Beta-1 protein (Recombinant) (Fc-Tag) (STJP000315)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered 100mM glycine PH 3.0.
Molecular Weight	Recombinant TGF-beta 1 is a protein consisting of 374 amino acid residue subunits. TGF-beta 1 1 has a predicted MW of 42kDa.
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C C for six months. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	7040
Gene Symbol	TGFB1
UniProt ID	TGFB1_HUMAN
Immunogen Region	Mature chain

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--