

Human IL2RA (CD25) protein (Recombinant) (STJP000332)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	0.2 Mu m filtered PBS solution, pH7.4.
Molecular Weight	Recombinant CD25 is a monomer protein consisting of 203 amino acid residue subunits, due to glycosylation migrates as an approximately 44 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	3559
Gene Symbol	IL2RA
UniProt ID	IL2RA_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



sIL-2R Alpha treated mice displayed increased novelty-induced stereotypic behaviors. In particular, sIL-2R Alpha treated mice showed dose-dependent increased head-bobbing