

Hepatitis B Virus surface antigen Adw2 protein (Recombinant) (STJP000337)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	P.pastoris

PRODUCT PROPERTIES

Formulation	20mM Na ₂ CO ₃ pH 9.2, 3mM EDTA, 15% glycerol.
Molecular Weight	Recombinant HBsAg migrates as an approximately 24 kDa protein under non-reducing conditions and reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Immunogen Region	Immunodominant Fragment
------------------	-------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--