

Human HIV-1 gp160 protein (Recombinant) (STJP000346)

GENERAL INFORMATION

Product Type	Proteins
--------------	----------

PRODUCT PROPERTIES

Formulation	0.2 μ m filtered 100mM PBS, pH4.5, 10mM tris HCL., 300mM Imidiasole, 8M Urea
Molecular Weight	Recombinant HIV gp160 is fused to an Beta-galactosidase purification tag and migrates as 58kDa on PAGE. NA
Storage Instruction	Can be stored in working aliquots at 4°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--