

Human CD66c (CEACAM6) protein (Recombinant) (STJP000369)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant Human CD66c is a monomer protein.The reduced monomer consists of 300 amino acids and has a predicted molecular mass of approx 32kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of CD66c monomer is approximately 50 k
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	1048
Gene Symbol	CEACAM5
UniProt ID	CEAM5_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--