

Mouse EPO protein (Recombinant) (STJP000374)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Mouse EPO is a protein consisting of amino acid residues [Ala 27-Arg192], due to glycosylation migrates as an approximately 27 kDa protein on SDS-PAGE. NA
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Gene ID	13856
Gene Symbol	Epo
UniProt ID	EPO_MOUSE
Immunogen Region	Mature chain
Immunogen Sequence	APPRLICDSR VLERYILEAK EAENVTMGCA EGPRLSENIT VPDTKVNFYA WKRMEEVEEQA IEVWQGLSL SEAILQAQAL LANSSQPPET LQLHIDKAIS GLRSLTSLR VLGAQKELMS PPDTPPPAPL RTLTVDTFCK LFRVYANFLR GKLLKLYTGEV NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--