

Human CSF-1 receptor protein (Recombinant) (STJP000375)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised 0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant human CSF-1 receptor is a monomer protein consisting of 504 amino acid residue subunits, due to glycosylation migrates as an approximately 90 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	1436
Gene Symbol	CSF1R
UniProt ID	CSF1R_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--