

Human LYN protein (Recombinant) (STJP000383)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Insect cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered 20mM Tris, 50mM NaCl, pH 8.0, 0.5mM PMSF.
Molecular Weight	Recombinant Human LYN is a monomer consisting of 736 amino acid residue, and migrates as an approximately 75 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	4067
Gene Symbol	LYN
UniProt ID	LYN_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--