

Mouse CTLA4 (CD152) protein (Recombinant) (STJP000404)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 M filtered PBS pH 7.4.
Molecular Weight	Recombinant Mouse CTLA4 (CD152) is a monomer consisting of 137 amino acid residue, and migrates as an approximately 25 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	12477
Gene Symbol	Ctla4
UniProt ID	CTLA4_MOUSE
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



Untransduced, immature DCs were treated on day 7 with either IFN-Gamma or CTLA4-Ig.