

## Mouse Lipocalin-2 protein (Recombinant) (STJP000411)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Host / Source | HEK293   |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.                                                                                                         |
| Molecular Weight    | Recombinant Mouse NGAL is a protein consisting of 185 amino acid residue subunits, due to glycosylation migrates as an approximately 25kDa protein on reduced SDS-PAGE. NA |
| Storage Instruction | Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA                                         |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">16819</a>      |
| Gene Symbol      | <a href="#">Lcn2</a>       |
| UniProt ID       | <a href="#">NGAL_MOUSE</a> |
| Immunogen Region | GLN20-ASP200, ECD          |

### ADDITIONAL INFORMATION

|      |                                                                                                              |
|------|--------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--------------------------------------------------------------------------------------------------------------|