

## Human Carboxypeptidase A1 protein (Recombinant) (STJP000433)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant Human Carboxypeptidase A1 (CPA1) is a protein consisting of 408 amino acid residue subunits, due to glycosylation migrates as an approximately as 50kDa band protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene ID	<a href="#">1357</a>
Gene Symbol	<a href="#">CPA1</a>
UniProt ID	<a href="#">CBPA1_HUMAN</a>
Immunogen Region	LYS17-TYR417

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--