

Human EPO receptor (EPOR) protein (Recombinant) (STJP000465)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant EPO receptor is a monomer protein consisting of 237 amino acid residue subunits, due to glycosylation migrates as an approximately 35 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	2057
Gene Symbol	EPOR
UniProt ID	EPOR_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--