

Human IL-32 protein (Recombinant) (STJP000485)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human IL32 is a protein consisting of amino acid residues [Met 1-Lys131], due to glycosylation migrates as an approximately 22 kDa protein on SDS-PAGE. NA
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Gene ID	9235
Gene Symbol	IL32
UniProt ID	IL32_HUMAN
Immunogen Region	ECD
Immunogen Sequence	MCFPKVLSDD MKKLKARMVM LLPTSAQGLG AWVSACDTE TVGHLGPWRD KDPALWCQLC LSSQHQAIER FYDKMQNAES GRGQVMSSLA ELEDDFKEGY LETVAAYYEE QHPELTPLE KERDGLRCRG NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--