

Mouse TNFR2 (CD120b) protein (Recombinant) (STJP000490)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	0.2 µg/ml filtered PBS solution, pH7.2.
Molecular Weight	Recombinant Mouse TNF receptor 2 (CD120b) is a protein consisting of 248 amino acid residue subunits, due to glycosylation migrates as an approximately 45kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	21938
Gene Symbol	Tnfrsf1b
UniProt ID	TNR1B_MOUSE
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--