

Human IL-3 receptor (CD123) protein (Recombinant) (STJP000494)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant IL-3RA (CD123) is a protein consisting of 298 amino acid residue subunits, due to glycosylation migrates as an approximately 50 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	3563
Gene Symbol	IL3RA
UniProt ID	IL3RA_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--