

Mouse IL-1 receptor type II (CD121b) protein (Recombinant) (STJP000500)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose, 0.01% tween80
Molecular Weight	Recombinant mouse IL-1 receptor type II (CD121b) is a protein consisting of 350 amino acid residue subunits, due to glycosylation migrates as an approximately 90 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	16178
Gene Symbol	Il1r2
UniProt ID	IL1R2_MOUSE
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--