

## Human CD30 (TNFRSF8) protein (Recombinant) (STJP000503)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Host / Source | HEK293   |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.                                                                                           |
| Molecular Weight    | Recombinant CD30 is a protein consisting of 370 amino acid residue subunits, due to glycosylation migrates as an approximately 80kDa protein on SDS-PAGE. NA |
| Storage Instruction | Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA                           |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">943</a>        |
| Gene Symbol      | <a href="#">TNFRSF8</a>    |
| UniProt ID       | <a href="#">TNR8_HUMAN</a> |
| Immunogen Region | ECD                        |

### ADDITIONAL INFORMATION

|      |                                                                                                              |
|------|--------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--------------------------------------------------------------------------------------------------------------|