

Mouse VEGFR2 protein (Recombinant) (Fc-Tag) (STJP000524)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Mouse KDR is a protein consisting of amino acid residues [Ala 20-Glu762], due to glycosylation migrates as an approximately 123 kDa protein on SDS-PAGE. NA
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Gene ID	16542
Gene Symbol	Kdr
UniProt ID	VGFR2_MOUSE
Immunogen Region	ECD
Immunogen Sequence	ASVGLPGDFL HPPKLSTQKD ILTILANTTL QITCRGQRDL DWLWPNAQRD SEERVLVTEC GGGDSIFCKT LTIPRVVGNL TGAYKCSYRD VDIASVYVY VRDYRSPFIA SVSDQHGIVY ITENKNKTVV IPCRGSISNL NVSLCARYPE KRFVPDGNRI SWDSEIGFTL PSYMISYAGM VFCEAKINDE TYQSIMYIVV VVGRIYDVI LSPPHEIELS AGEKLVLN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--