

Human TGF Beta-1 receptor 2 (TGRBR2) protein (Recombinant) (STJP000530)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered 100mM Glycine, 10mM NaCl, 50mM Tris, pH 7.5.
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C C for six months. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	7048
Gene Symbol	TGFB2
UniProt ID	TGFR2_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--