

Mouse 4-1BB Ligand (TNFSF9) protein (Recombinant) (STJP000540)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised 0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant Mouse 4-1BBL is a protein consisting of 220 amino acid residue subunits, due to glycosylation migrates as an approximately 25kDa protein on SDS-PAGE under reducing conditions. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	21950
Gene Symbol	Tnfsf9
UniProt ID	TNFL9_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--