

## Human IGF-1R (CD221) protein (Recombinant) (STJP000574)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant human IGF-1R (CD221) is a heterodimer protein, due to glycosylation migrates as an approximately 50kDa, 120kDa, 150kDa protein under reducing ocnditions on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year, addition of a carrier protein to your working aliquots (0.1% BSA or HSA) is recommended for long term storage Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Immunogen Region	ECD
------------------	-----

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--