

Mouse CD80 (B7-1) protein (Recombinant) (STJP000583)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2um filtered 50mM Tris, 100mM NaCl, pH 8.0.
Molecular Weight	Recombinant mouse CD80 is a monomer protein consisting of 220 amino acid residue subunits, due to glycosylation migrates as an approximately 45 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	12519
Gene Symbol	Cd80
UniProt ID	CD80_MOUSE
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--