

Human BCMA (CD269) protein (Recombinant) (STJP000587)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M unfiltered mM Tris solution pH7.5.
Molecular Weight	Recombinant BCMA (CD269) Fc chimera is a homodimer protein consisting of 300 amino acid residue subunits, and migrates as an approximately 80kDa band under non-reducing and as a 40Da band under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 4°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	608
Gene Symbol	TNFRSF17
UniProt ID	TNR17_HUMAN
Immunogen Region	ECD
Immunogen Sequence	MLQMAGQCSQ NEYFDSLHA CIPCQLRCSS NTPPLTCQRY CNASVTNSVK GTNA ILWTCL NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--