

Mouse PD-1 (PDCD1) protein (Recombinant) (STJP000588)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant Mouse PD-1 (PDCD1) is a protein consisting of 147 amino acid residue subunits, due to glycosylation migrates as an approximately 40kDa protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	18566
Gene Symbol	Pdc1
UniProt ID	PDCD1_MOUSE
Immunogen Region	LEU25-GLN167

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--