

## Rat IL-1 receptor 1 (CD121a) protein (Recombinant) (STJP000591)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant IL-R1 (CD121a) is a protein consisting of 330 amino acid residue subunits, due to glycosylation migrates as an approximately 55-60 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene Symbol	<a href="#">Il1r1</a> Uniprot ID= <a href="#">IL1R1_RAT</a> "
Immunogen Region	ECD

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--