

Mouse CD31 (PECAM1) protein (Recombinant) (STJP000599)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant mouse CD31 (PECAM-1) is a monomer protein consisting of 585 amino acid residue subunits, due to glycosylation migrates as an approximately 100 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	5175
Gene Symbol	PECAM1
UniProt ID	PECA1_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--