

Human CD157 (BST1) protein (Recombinant) (STJP000600)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	0.2 µg/ml filtered PBS solution, pH7.2.
Molecular Weight	Recombinant CD157 is a monomer protein consisting of 275 amino acid residue subunits, due to glycosylation migrates as an approximately 38-43 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	683
Gene Symbol	BST1
UniProt ID	BST1_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--