

## Human ADAM15 protein (Recombinant) (STJP000603)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant human ADAM15 (MDC15) is a protein consisting of 497 amino acid residue subunits, due to glycosylation migrates as an approximately as 70kDa band protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene ID	<a href="#">8751</a>
Gene Symbol	<a href="#">ADAM15</a>
UniProt ID	<a href="#">ADA15_HUMAN</a>
Immunogen Region	ASP207-THR696

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--------------------------------------------------------------------------------------------------------------