

Human IL-35 protein (Recombinant) (STJP000607)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO Cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution pH 7.4.
Molecular Weight	Recombinant Human Interleukin-35 is a monomer consisting of 645 amino acid residue, and migrates as an approximately 95 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	10148
Gene Symbol	EBI3
UniProt ID	IL27B_HUMAN
Immunogen Region	Mature chain

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--