

Human IL-13 alpha 2 receptor (IL13RA2) protein (Recombinant) (STJP000610)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised 0.2 Mu m filtered 50mM Tris, 100mM glycine solution, pH7.5 and 5% trehalose.
Molecular Weight	Recombinant IL-13RA2 is a homodimer protein consisting of 563 amino acid residue subunits, due to glycosylation migrates as an approximately 130kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	3598
Gene Symbol	IL13RA2
UniProt ID	I13R2_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--