

Hemagglutinin Influenza A-Indonesia-05-2005-H5N1 protein (Recombinant) (STJP000616)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Insect cells

PRODUCT PROPERTIES

Formulation	Lyophilised from PBS pH7.4, 5% trehalose, 5% mannitol.
Molecular Weight	Recombinant InfluenzaH5N1 Hemagglutinin is a monomer protein, and migrates due to glycosylation as an approximately 70 kDa protein under non-reducing conditions and as 45kDa and 25kDa fragments under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Immunogen Region	ASP17-GLN531
------------------	--------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--