

## Mouse IL-17RA (CD217) protein (Recombinant) (STJP000625)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised 0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant Human IL18RA is a monomer protein consisting of 322 amino acid residue subunits, due to glycosylation migrates as an approximately 65kD protein under reducing SDS-PAGE conditions. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene ID	<a href="#">16172</a>
Gene Symbol	<a href="#">Il17ra</a>
UniProt ID	<a href="#">I17RA_MOUSE</a>
Immunogen Region	ECD

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--