

## Mouse CD33 (Siglec-3) protein (Recombinant) (STJP000633)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant CD33 (Siglec-3) is a Homodimer protein consisting of two 460 amino acid residue subunits, due to glycosylation migrates as an approximately 60kDa protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene ID	<a href="#">12489</a>
Gene Symbol	<a href="#">Cd33</a>
UniProt ID	<a href="#">CD33_MOUSE</a>
Immunogen Region	ECD

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--