

Mouse IL-11 receptor alpha protein (Recombinant) (STJP000637)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution pH7.4.
Molecular Weight	Recombinant Mouse Interleukin-11RA is a monomer consisting of amino acid residue subunits, and migrates as an approximately 46 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 4°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	16157
Gene Symbol	Il11ra1
UniProt ID	I11RA_MOUSE
Immunogen Region	ECD
Immunogen Sequence	MSSSCSGLTR VLVAVATALV SSSSPCPQAW GPPGVQYGQP GRPVMGCCPG VSAGTPVSWF RDGDSRLQGPDSGLGHRV LAQVDSPDEG TYVCQTLQGV SGGMVTLLKGFPPARPEVSC QAVDYENFSC TWSPGQVSGLPTRYLTSYRK KTLPGAESQR ESPSTGPWPC PQDPLEASRC VHGAEFWSE YRINVTEVNP LGASTCLLDV RLQSILRPDP PQGLRVES

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--