

## Human IL-20R2 protein (Recombinant) (Fc-Tag) (STJP000642)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised 0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant Human IL-20RB is a homodimer protein consisting of 210 amino acid residue subunits, due to glycosylation migrates as an approximately 100kD protein under reducing SDS-PAGE conditions. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene ID	<a href="#">53833</a>
Gene Symbol	<a href="#">IL20RB</a>
UniProt ID	<a href="#">I20RB_HUMAN</a>
Immunogen Region	ASP30-ALA230

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--