

Mouse Carnosinase protein (Recombinant) (STJP000649)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant Mouse Carnosinase (CN1) is a protein consisting of 500 amino acid residue subunits, due to glycosylation migrates as an approximately as 55kDa band protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	338403
Gene Symbol	Cndp1
UniProt ID	CNDP1_MOUSE
Immunogen Region	MET1-TYR492

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--