

Mouse CD155 protein (Recombinant) (STJP000652)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	0.2 µg/ml filtered PBS solution, pH7.2.
Molecular Weight	Recombinant mouse PVR (CD155) is a monomer protein consisting of 329 amino acid residue subunits, due to glycosylation migrates as an approximately 55-60 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Immunogen Region	ECD
------------------	-----

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--