

Human IL-17RC protein (Recombinant) (Fc-Tag) (STJP000656)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised 0.2 Mu m filtered PBS solution, pH7.2.
Molecular Weight	Recombinant IL17RC-Fc is a homodimer protein consisting of 675 amino acid residue subunits, due to glycosylation migrates as an approximately 100kD protein under reducing SDS-PAGE conditions. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	84818
Gene Symbol	IL17RC
UniProt ID	I17RC_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--