

Rat Interleukin-12 p40 protein (Recombinant) (STJP000666)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Insect cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution.
Molecular Weight	Native rat Interleukin-12 is generated by the proteolytic removal of the signal peptide and propeptide, the molecule has a calculated molecular mass of approximately 70 kDa. Recombinant rat Interleukin-12 is a disulfide-linked heterodimeric protein,
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Immunogen Region	Mature chain
Immunogen Sequence	NA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--